

Anti-Serum Amyloid P/SAP 兔多克隆抗体

K1337863 -20°C



产品信息

蛋白质全称	血清淀粉样蛋白 P 组分
别名	SAP, Apcs, Ptx2, Amyloid P component serum
免疫原	KLH 偶联的小鼠 Serum Amyloid P/SAP 合成肽
Uniprot ID	P12246
抗体亚型	IgG
纯化方式	亲和纯化
预测分子量 / 观测分子量	26 kDa / 24 kDa

应用

WB 免疫印迹

物种：人, 小鼠
 稀释比：1: 500-1: 1000
 阳性样品：293T 细胞


背景

The protein encoded by this gene is a glycoprotein, belonging to the pentraxin family of proteins, which has a characteristic pentameric organization. These family members have considerable sequence homology which is thought to be the result of gene duplication. The binding of the encoded protein to proteins in the pathological amyloid cross-beta fold suggests its possible role as a chaperone. This protein is also thought to control the degradation of chromatin. It has been demonstrated that this protein binds to apoptotic cells at an early stage, which raises the possibility that it is involved in dealing with apoptotic cells in vivo.

储存条件

储存条件：自收到之日起 -20°C下储存一年，避免反复冻融。
储存缓冲液：含 0.15% ProClin300 防腐剂, 100 μg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



WB 检测 Serum Amyloid P/SAP 蛋白(货号 K1337863).
 样品：经 RIPA 裂解液(货号 KSG2002)处理的蛋白质。
 封闭：3%脱脂牛奶(货号 KSGC310001)溶于 TBST 溶液，室温孵育 1 小时。
 一抗：1: 1000 稀释, 4°C 孵育过夜。
 二抗：HRP 标记山羊抗兔 IgG (H+L) (货号 KB69909), 1: 5000 稀释, 室温孵育 1 小时。

For Research Use Only!

Ver. No.: V1-202603