

Anti-Septin 8 兔多克隆抗体

K1338319 -20°C



产品信息

蛋白质全称	Septin-8
别名	Septin8, Kiaa0202, Sept8, SEP2
免疫原	小鼠 Septin 8 重组蛋白
Uniprot ID	Q8CHH9
抗体亚型	IgG
纯化方式	亲和纯化
预测分子量 / 观测分子量	50 kDa / 50 kDa

应用

WB 免疫印迹

物种：人, 小鼠, 大鼠

稀释比：1: 3000-1: 6000

阳性样品：心

背景

SEPT8, also named as KIAA0202, is a member of the highly conserved septin family. Septins are 40-60kD GTPases that assemble as filamentous scaffolds. They are involved in the organization of submembranous structures, in neuronal polarity, and in vesicle trafficking. It may play a role in platelet secretion. SEPT8 is about 50-58 kDa in WB detection.

储存条件

储存条件：自收到之日起 -20°C下储存一年，避免反复冻融。
储存缓冲液：含 0.15% ProClin300 防腐剂, 100 μg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



WB 检测 Septin 8 蛋白(货号 K1338319).
 样品：经 RIPA 裂解液(货号 KSG2002)处理的蛋白质。
 封闭：3%脱脂牛奶(货号 KSGC310001)溶于 TBST 溶液，
 室温孵育 1 小时。
 一抗：1: 3000 稀释, 4°C 孵育过夜。
 二抗：HRP 标记山羊抗兔 IgG (H+L) (货号 KB69909), 1:
 5000 稀释, 室温孵育 1 小时。

For Research Use Only!

Ver. No.: V1-202603