

Anti-PAX7 兔多克隆抗体

K1339570 -20°C



产品信息

蛋白质全称	配对盒蛋白 Pax7
别名	HUP1, paired box 7, Paired box protein Pax 7, PAX7, PAX7B, RMS2, PAX7 transcriptional factor
免疫原	小鼠 PAX7 重组蛋白
Uniprot ID	P47239
抗体亚型	IgG
纯化方式	亲和纯化
亚细胞定位	细胞核

应用

IHC/IF 免疫组织化学/免疫荧光

物种：小鼠

稀释比：1: 2100-1: 4200

阳性样品：胎脑


背景

PAX7, also named as HUP1, belongs to the paired homeobox family. It is probable a transcription factor. PAX7 may have a role in myogenesis. A chromosomal aberration involving PAX7 is a cause of rhabdomyosarcoma 2 (RMS2) which also known as alveolar rhabdomyosarcoma. The antibody is specific to PAX7. It will not cross react with other PAX members.


储存条件

储存条件：自收到之日起 -20°C下储存一年，避免反复冻融。
储存缓冲液：含 0.15% ProClin300 防腐剂, 100 µg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



IHC 检测 PAX7 蛋白(货号 K1339570).
 样品: 小鼠胎脑, 4%多聚甲醛 (货号 KSG1101) 固定 12-24 小时.
 抗原修复: 柠檬酸抗原修复液(干粉, pH 6.0) (KSG1201), 98°C, 20 分钟.
 一抗: 1: 3200 稀释, 4°C 孵育过夜.
 二抗: S-vision 免疫组化多聚二抗(山羊抗兔), 即用型 (货号 KB3906), 室温孵育 20 分钟.



IHC 检测 PAX7 蛋白(货号 K1339570).
 样品: 小鼠胎脑, 4%多聚甲醛 (货号 KSG1101) 固定 12-24 小时.
 抗原修复: 柠檬酸抗原修复液(干粉, pH 6.0) (KSG1201), 98°C, 20 分钟.
 一抗: 1: 3200 稀释, 4°C 孵育过夜.
 二抗: S-vision 免疫组化多聚二抗(山羊抗兔), 即用型 (货号 KB3906), 室温孵育 20 分钟.

For Research Use Only!

Ver. No.: V1-202603