

Anti-TEX11 兔多克隆抗体

K1340062 -20°C



产品信息

蛋白质全称	睾丸表达蛋白 11
别名	testis expressed 11, TEX11, TGC1, TSGA3, ZIP4H, Zip4, Protein ZIP4 homolog
免疫原	小鼠 TEX11 重组蛋白
Uniprot ID	Q14AT2
抗体亚型	IgG
纯化方式	亲和纯化
亚细胞定位	细胞核

应用

IHC/IF 免疫组织化学/免疫荧光

物种: 小鼠, 大鼠

稀释比: 1: 1000-1: 2000

阳性样品: 睾丸


背景

Testis-expressed 11 (TEX11) is an X-chromosome-encoded germ-cell-specific protein that is expressed most abundantly in spermatogonia and early spermatocytes in the testes. The TEX11 protein was found to form discrete foci along the synaptonemal complexes (SCs) on meiotic chromosomes and was shown to be essential for meiosis in males.


储存条件

储存条件: 自收到之日起 -20°C下储存一年, 避免反复冻融。
储存缓冲液: 含 0.15% ProClin300 防腐剂, 100 µg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



IHC 检测 TEX11 蛋白(货号 K1340062).
 样品: 小鼠睾丸, 4%多聚甲醛 (货号 KSG1101) 固定 12-24 小时.
 抗原修复: Tris-EDTA 抗原修复液(pH 9.0) (KSG1203), 98°C, 20 分钟.
 一抗: 1: 1500 稀释, 4°C 孵育过夜.
 二抗: S-vision 免疫组化多聚二抗(山羊抗兔), 即用型 (货号 KB3906), 室温孵育 20 分钟.



IHC 检测 TEX11 蛋白(货号 K1340062).
 样品: 大鼠睾丸, 4%多聚甲醛 (货号 KSG1101) 固定 12-24 小时.
 抗原修复: Tris-EDTA 抗原修复液(pH 9.0) (KSG1203), 98°C, 20 分钟.
 一抗: 1: 1500 稀释, 4°C 孵育过夜.
 二抗: S-vision 免疫组化多聚二抗(山羊抗兔), 即用型 (货号 KB3906), 室温孵育 20 分钟.

For Research Use Only!

Ver. No.: V1-202603