

Anti-REA 兔多克隆抗体

K1341178 -20°C



产品信息

蛋白质全称	抗增殖蛋白 2
别名	BAP, Bap37, D prohibitin, p22, PHB2, PNAS 141, REA, Repressor of estrogen receptor activity, B-cell receptor-associated protein BAP37
免疫原	KLH 偶联的小鼠 REA 合成肽
Uniprot ID	O35129
抗体亚型	IgG
纯化方式	亲和纯化
预测分子量 / 观测分子量	33 kDa / 33 kDa

应用

WB 免疫印迹

物种：人, 小鼠, 大鼠

稀释比：1: 700-1: 1400

阳性样品：心


背景

PHB2, also named as BAP, REA, D-prohibitin, BAP37 and Prohibitin-2, belongs to the prohibitin family. It acts as a mediator of transcriptional repression by nuclear hormone receptors via recruitment of histone deacetylases. PHB2 functions as an estrogen receptor (ER)-selective coregulator that potentiates the inhibitory activities of antiestrogens and represses the activity of estrogens. It competes with NCOA1 for modulation of ER transcriptional activity. PHB2 is probably involved in regulating mitochondrial respiration activity and in aging.

储存条件

储存条件：自收到之日起 -20°C下储存一年，避免反复冻融。
储存缓冲液：含 0.15% ProClin300 防腐剂, 100 µg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



WB 检测 REA 蛋白(货号 K1341178).
 样品: 经 RIPA 裂解液(货号 KSG2002)处理的蛋白质.
 封闭: 3%脱脂牛奶(货号 KSGC310001)溶于 TBST 溶液,
 室温孵育 1 小时.
 一抗: 1: 1000 稀释, 4°C 孵育过夜.
 二抗: HRP 标记山羊抗兔 IgG (H+L) (货号 KB69909), 1:
 5000 稀释, 室温孵育 1 小时.

For Research Use Only!

Ver. No.: V1-202603