

Anti-DUSP4 兔多克隆抗体

K1335733 -20°C



产品信息

蛋白质全称	双特异性蛋白磷酸酶 4
别名	dual specificity phosphatase 4, DUSP4, HVH2, MAP kinase phosphatase 2, MKP 2, MKP2, MKP-2, TYP, VH2
免疫原	小鼠 DUSP4 重组蛋白
Uniprot ID	Q8BFV3
抗体亚型	IgG
纯化方式	亲和纯化
亚细胞定位	细胞核

应用

IHC/IF 免疫组织化学/免疫荧光

物种：人
 稀释比：1: 1000-1: 2000
 阳性样品：脑


背景

DUSP4, also named as MKP-2, is a member of the mitogen-activated protein (MAP) kinase phosphatase family which has been suggested to play an important role in the feedback control of MAP kinase-mediated gene expression. DUSP4 preferentially inactivates ERK1, ERK2 and JNK, is expressed in a variety of tissues, and is localized in the nucleus. DUSP4 encodes two distinct isoforms: MKP-2-L(43 kDa) and MKP-2-S (33 kDa).


储存条件

储存条件：自收到之日起 -20°C下储存一年，避免反复冻融。
储存缓冲液：含 0.15% ProClin300 防腐剂, 100 μg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



IHC 检测 DUSP4 蛋白(货号 K1335733).
 样品: 人脑, 4%多聚甲醛 (货号 KSG1101) 固定 12-24 小时.
 抗原修复: 柠檬酸抗原修复液(干粉, pH 6.0) (KSG1201), 98°C, 20 分钟.
 一抗: 1: 1500 稀释, 4°C 孵育过夜.
 二抗: S-vision 免疫组化多聚二抗(山羊抗兔), 即用型 (货号 KB3906), 室温孵育 20 分钟.



IHC 检测 DUSP4 蛋白(货号 K1335733).
 样品: 人脑, 4%多聚甲醛 (货号 KSG1101) 固定 12-24 小时.
 抗原修复: 柠檬酸抗原修复液(干粉, pH 6.0) (KSG1201), 98°C, 20 分钟.
 一抗: 1: 1500 稀释, 4°C 孵育过夜.
 二抗: S-vision 免疫组化多聚二抗(山羊抗兔), 即用型 (货号 KB3906), 室温孵育 20 分钟.

For Research Use Only!

Ver. No.: V1-202603