

Anti-TSPAN15 兔多克隆抗体

K1341667 -20°C



产品信息

蛋白质全称	四次跨膜蛋白 15
别名	NET-7, NET7, Tetraspan NET 7, tetraspanin 15, TM4SF15, Tspan-15, TSPAN15, 2700063A19Rik
免疫原	小鼠 TSPAN15 重组蛋白
Uniprot ID	F7BWT7
抗体亚型	IgG
纯化方式	亲和纯化
预测分子量 / 观测分子量	33 kDa / 33 kDa

应用

WB 免疫印迹

物种：人, 小鼠, 大鼠

稀释比：1: 1000-1: 2000

阳性样品：胸腺

背景

The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The use of alternate polyadenylation sites has been found for this gene.

储存条件

储存条件：自收到之日起 -20°C下储存一年，避免反复冻融。
储存缓冲液：含 0.15% ProClin300 防腐剂, 100 μg/mL 牛血清白蛋白和 50%甘油的磷酸盐缓冲液。

图像



WB 检测 TSPAN15 蛋白(货号 K1341667).
 样品：经 RIPA 裂解液(货号 KSG2002)处理的蛋白质。
 封闭：3%脱脂牛奶(货号 KSGC310001)溶于 TBST 溶液，室温孵育 1 小时。
 一抗：1: 1000 稀释, 4°C 孵育过夜。
 二抗：HRP 标记山羊抗兔 IgG (H+L) (货号 KB69909), 1: 5000 稀释, 室温孵育 1 小时。

For Research Use Only!

Ver. No.: V1-202603